# E 315/1
## Full demineralisation single use resin
## For the regeneration of full demineralisation cartridges

EAN: 4002516957515 / Material number: 12945010 / Old Material Number: 69431502D



- For use in the LP 2800/1 empty cartridge
- 20 l of homogeneously mixed resins
- Carton with 2 x 10 l bags
- Vacuum sealed in plastic bags
- Including replacement filter bag

| **Product category** | |
|---|---|
| Fresh-water dishwashers | • |
| Tank dishwashers | • |
| Throughfeed tank dishwashers | • |
| Laboratory washers | • |
| Thermal disinfectors | • |
| Washer-disinfectors, medical | • |
| **Machine affiliation** | |
| PG 8067 | • |
| PG 8041 | • |
| PG 8043 | • |
| PG 8045 | • |
| PG 8055 | • |
| PG 8055 U | • |
| PG 8056 | • |
| PG 8056 U | • |
| PG 8057 TD | • |
| PG 8057 TD U | • |
| PG 8058 | • |
| PG 8058 U | • |
| PG 8059 | • |
| PG 8059 U | • |
| PG 8060 | • |
| PG 8061 | • |
| PG 8061 U | • |
| PG 8063 | • |
| PG 8164 | • |
| PG 8165 | • |
| PG 8166 | • |
| PG 8169 | • |
| PG 8172 | • |
| PG 8172 Eco | • |
| PTD 901 | • |
| PG 8096 | • |
| PG 8096 U | • |
| PG 8099 | • |
| PG 8099 U | • |
| PFD 400 | • |
| PFD 400 U | • |
| PFD 401 | • |
| PFD 401 U | • |
| PFD 402 | • |
| PFD 402 U | • |
| PFD 404 | • |
| PFD 404 U | • |
| PFD 405 | • |
| PFD 405 U | • |
| PFD 407 | • |

E 315/1
Full demineralisation single use resin
For the regeneration of full demineralisation cartridges

EAN: 4002516957515 / Material number: 12945010 / Old Material Number: 69431502D



| | |
|---|---|
| PFD 407 U | • |
| PFD 408 U | • |
| PTD 701 | • |
| PTD 702 | • |
| PTD 703 | • |
| PTD 704 | • |
| PLW 8683 | • |
| PLW 8683 CD | • |
| PWD 8682 | • |
| PWD 8682 C | • |
| PWD 8692 | • |
| **Product category** | |
| Water treatment | • |
| **Dimensions and weight** | |
| External dimensions, gross height in mm | 285 |
| External dimensions, gross width in mm | 260 |
| External dimensions, gross depth in mm | 440 |
| Net weight in kg | 13.77 |
| Gross weight in kg | 14.27 |